

FAQ – NovaPro™ Flow

1. What is NovaPro™ Flow?

A nanofiber reinforced, nano-hybrid flowable composite for the use of dental restorations.

2. What quantities and unit sizes are offered?

Available in 1 and 2 syringe packs with 10 disposable tips/syringe. Each syringe contains 2.0 grams of composite. We also offer an Intro Kit of our most popular shades (A1, A2, A3.5, and B1).

3. What shade options are available for NovaPro™ Flow?

Available in A1, A2, A3, A3.5, B1, B2, B3, and C1 shades. All our shades follow the VITA Classical Shade Guide.

4. Does NovaPro™ Flow contain any allergens or toxins?

No major allergens or toxins (i.e., mercury, nuts, gluten, shellfish, dairy, wheat, and soy).

5. What is the shelf life?

2 years. We recommend storing in a dry, cool location for maximum shelf life (i.e., refrigerator).

6. Which types/classes can NovaPro™ Flow restore?

Class I, II, III, and V. It can also be used as a base/liner or sealant.

7. What is the cure time and depth recommendations?

10 second cure time. If you prefer the standard 20 seconds, there will be no negative effects. The composite should be applied in 2mm increments.

8. Are there any differences in using our flowable compared to other flowables?

No. This product can be used the same way you use any other flowable composite.

9. Why switch to NovaPro™ Flow? How many dentists are making the switch?

We are the only composite that contains patented nanofiber technology, allowing for longer lasting restorations, higher flexural strength, and high percentage of cure, leading to less re-care and breakdown of restorations. Many dentists are choosing NovaPro™ based on these benefits.

10. Where is NovaPro™ Flow made and sold?

We proudly manufacture and sell all of our products nationwide from our facility in Columbia, Missouri. Our nanofiber technology was developed in collaboration with the University of Missouri.